

Product Name :
Heating Oven**Product Code :**
NMCBS11087**Description :**

We offer laboratory ovens that are used in a variety of thermal processing applications including general lab work, component and stability testing, core hardening, drying glassware, and sterilizing. These are known for their robust construction, low maintenance and high efficiency. Some of the features are as follows: •Double Walled Construction •Exterior & Interior made of stainless Heating Oven steel 304/316 (GMP model) or Exterior made of Mild Steel with powder coating finish and interior made of stainless steel 304 (STD Model) •80/20 Ni Heating Elements placed at bottom •Temp. Range : 5°C above ambient to 300°C 2°C Heating Oven, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,