

**Product Name :**  
Curing Tank**Product Code :**  
NMCBS11090**Description :**

Curing Tank We offer Curing Tanks that are used in a variety of thermal processing applications including general laboratory work, component and stability testing, drying glassware, and sterilizing. These are available in various specifications thus meeting the varied needs of our clients. Some of the features are as follows: •Double Walled – Outer MS – Inner MS or St.Steel. •80/20 Ni Heating Elements placed at two vertical sides in addition to the bottom one for uniform heating in the chamber. •Fan arrangement for forced air circulation through out the chamber. •Temp. Range : 5°C above ambient to 300°C 2°C •230V AC 50HZ Mains Supply. Curing Tank, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,