

**Product Name :**  
Accelerated Curing Tank**Product Code :**  
NMCBS11091**Description :**

Accelerated Curing Tank We are manufacturers of Accelerated Curing ovens that are used in a variety of thermal processing applications including general lab work, component and stability testing, core hardening and sterilizing. These are available in various specifications thus meeting the varied needs of our clients. Some of the features are as follows: •Double Walled – Outer MS – Inner MS or St.Steel. •80/20 Ni Heating Elements placed at two vertical sides in addition to the bottom one for uniform heating in the chamber. •Temp. Range : 5°C above ambient to 300°C 2°C •230V AC 50HZ Mains Supply. •Digital Temperature Controller. Accelerated Curing Tank, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers