

**Product Name :**  
Seed Germinator**Product Code :**  
NMCBS11095**Description :**

Seed Germinator Double walled-Inner SS/Outer MS 75 mm Gap filled with puff insulation; Digital display preset electronic controllers for temperature and dew point. Hermetically sealed compressor unit for cooling system; castor wheels for easy movability Pilot lamps and switches for all controls 230V AC 50HZ Mains. Temp. Range: 10°C to 60°C±1°C, Humidity: 40% to 95% Rh±3% Rh. Direct Rh Indicator SS water tank is provided with a ball cock valve & front door is provided with a full view glass window for inner observation. Seed Germinator, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,