



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Semi Fowler Hospital Bed

Product Code :
NMHOPB00003



Description :

Semi Fowler Hospital Bed Size : 198 L x 90 W cm x 60 H cms.

Technical Specification :

Semi fowler bed made of SS
Suitable to be used in hospitals and Covid centers to treat patients by giving them head rest
Easy to move with elevation angle from 0 to 85 degree

Note : Semi Fowler Hospital Bed can also supplied with castors.



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,