

**Product Name :**  
U.V. Cabinet**Product Code :**  
NMCBS11096**Description :**

U.V. Cabinet solutions offered find usage for conducting chromatographic analysis through use of ultra violet fluorescence under laboratory lightning conditions having independent control at the top which makes the operation fast simple as well as convenient to handle. As a self contained unit, the cabinet is designed to view as well as make chromatograms and be also used in general fluorescence studies. Further, the cabinet contains long wave, short wave ultra violet light as well as white light lamp in compact table top that helps in providing for three means of locating as well as identifying of components. These chromatograms can be viewed with 3600-4000 A long waves or as 2500 A short waves ultra violet light or white light. Here, the presence of black curtain ensures darkness as well as brilliant lightning in cabinet and at same time also allows users to handle sample and switches with complete ease. For durable usage, the cabinet is constructed using sturdy metal with two side handles so as to provide for easy transportation. Further, the switches are also fitted on top for easy of handling while handling detection of traces of minerals and studying of old inspection U.V. Cabinet, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](#)