

Product Name :
Vertical Autoclave**Product Code :**
NMCBS11097**Description :**

Vertical Autoclave solutions offered find application in day-to-day sterilization process in different laboratory applications and feature use of both standard and customized models that are scientifically designed so as to meet demands of scientists, individuals as well as for specialized research applications. Here, our expertise also lies in offering economy models to many reputed customers in India including those in Defence sector, for Research laboratories, Educational institutes, R & D labs as well as big national and international companies. For achieving durable performance, these vertical autoclaves feature double walled construction with the outer wall made of thick PCRC sheet that are pre-treated with primer and rust proofing done with lasting stove enamel or through powder coating. Other than this, heavy gauge stainless steel of SS 304 grade comprises inner wall with mineral glass wool insulation provided on customers' demand. Further, these systems deliver heating through specially designed water heaters of defined wattage so as to ensure optimum sterilization time in every sterilization cycle. These Vertical Autoclaves have working temperature of 121°C as per defined international standards and guidelines which can be increased to 134°C based on Vertical Autoclave, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,