

Product Name :
Refrigerated Water Bath**Product Code :**
NMCBS11109**Description :**

Refrigerated Water Bath Refrigerated Water Bath is used for carrying out day to day bacteriological as well as varied laboratory applications that require incubation and general tests and procedures including Wassermann, Kilmer tests as followed in research and testing laboratories. Further, refrigerated water baths also have usage in field of carrying out enzyme reaction studies; tissues culture application, doing growth observation studies, arriving at fermentation analysis as well as handling other general specialized applications in a laboratory environment. Here, our expertise lies in offering these in model options that meet the challenging demands of scientists and of people engaged in conducting specialized research. For durable performance, these are also fitted with ISI make compressor so as to ensure effective cooling and working temperature range from 5°C to 90°C. Applications: Used for Beer Cards tank test. Accord with standard "GB/T4928-2001". Top market share in China market for 5 years. Salient Features: Ergonomic Design. Energy Efficient. Long life. Low maintenance. Calibration and protocol documentation. Refrigerated Water Bath, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers