

Product Name :
Mortuary Freezer**Product Code :**
NMCBS11115**Description :**

Mortuary Freezer Owing to our vast experience, we have been successful in catering to the requirements of our esteemed clients by offering qualitative range of Mortuary Freezer. These products are used for storing the cadaverous under low temperature cool condition in order to prevent decomposition. These products are double walled with the outer surface made of M.S. sheet which is duly painted and inner wall made of stainless steel. Mortuary Freezer front opening covered door is provided with the magnetic gasket lock & key and handles individual chambers provided with independent lock and duplicate keys. Technical Details: · Made of G.I. sheet so as to carry dead bodies and are provided with handles · Trays travel on rollers fitted with ball bearing for easy movement · Hermetically sealed compressor serves refrigeration unit for maintaining the temperature range of 2°C to 5°C or (-) 15°C to (-) 20°C (as ordered) · Unit is fitted with Digital Temperature Indicator-cum-controller · Works on 220/230V AC. Power supply Specification: Capacity Temp. Temp. Voltage Stabilizer Single Body 2° to 5° C -15° to -20° C 3 KVA Two Bodies 2° to 5° C -15° to -20° C 4 KVA Three Bodies 2° to 5° C -15° to -20° C 5 KVA Four Bodies 2° to 5° C -15° to -20° C 8 KVA Six Bodies 2° to 5° C -15° to -20° C 8 KVA Mortuary Freezer, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers