

**Product Name :**  
Precision Water Bath**Product Code :**  
NMCBS11119**Description :**

Precision Water Bath Keeping in sync with the times, we also offer a wide range of Precision Water Bath. These are also called as constant temp. control. Our products are made as per ADI-57 and are provided with glass windows on both opposite sides for observation. Precision Water Bath temperature of these products are controlled by mercury contact thermometer with a sensitivity of  $\pm 1^{\circ}\text{C}$ . Our products are suitable for precision laboratory activities requiring constant temperature with high degree of accuracy. These are a complete stirrer and are workable on 220V AC 50 Hz single phase. Optional: · Digital temperature indicator-cum-controller instead of contact thermometer · Microprocessor PID digital temp. indicator-cum-controller instead of contact thermometer Features: · Sturdy construction · Compact design · High performance Specifications: Chamber size in mm & inches. 425x250x225 mm (17"x10"x9") 450x375x375 mm (18"x15"x15") 600x450x450 mm (24"x18"x18") Precision Water Bath, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,