

**Product Name :**  
Solid Liquid Extraction

**Product Code :**  
NMCBS11138



**Description :**

Solid Liquid Extraction With firm commitment to quality, we are engaged in offering excellent quality, Solid Liquid Extraction (Packed Bed Type) . These are designed under the strict command of our quality analyzers keeping in mind the industry guidelines. Post production, these products are stringently tested on varied grounds, to ensure their compliance with international quality standards. Features: Durability High tensile strength Resistance to corrosion Solid Liquid Extraction, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers