

**Product Name :**  
Pelton Wheel Turbine Test Rig**Product Code :**  
NMCBS11168**Description :**

Pelton Wheel Turbine Test Rig Our clients can avail from us high quality Pelton Wheel Turbine Test Rig to the esteemed customer, as per their specific requirements. Widely used for experimental purposes, these products are highly demanded in various laboratories and industries. The range is manufactured using cast iron, steel or brass. Moreover, these products are fitted with a big transparent acrylic window, nozzle & spear arrangement, rotor assembly of shafts, runner and brake drum that are mounted on a strong base frame. Features: Available at a very exciting and reasonable range of prices Very efficient in its performance Long lasting and durable Robust structure Hassle free functioning Easy to operate Accurate value of reading Pelton Wheel Turbine Test Rig, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](#)