

**Product Name :**  
Semi Automatic Rotary Microtome**Product Code :**  
NMCBS13010**Description :**

Semi Automatic Rotary Microtome We are offering Semi Automatic Rotary Microtome catering to the needs and requirements of our esteemed customers. This equipment is featured with a standard drive system for smooth sectioning and trimming. We have also installed LCD screen displays section, various trimming thickness range and number of sections in this microtome. There is a lockable hand wheel, interchangeable specimen clamp along with a safety alarm system and huge volume waste tray for easy operation. Semi Automatic Rotary Microtome calibrated controls assure an accurate positioning to make re-cuts possible. Specifications: · Section thickness setting range: 0.25-60  $\mu\text{m}$  · Trimming thickness setting Range: 0.25-60  $\mu\text{m}$  · Setting values: 0.25-1  $\mu\text{m}$  increment step 0.25  $\mu\text{m}$  · 1-10  $\mu\text{m}$  increment step 1  $\mu\text{m}$  · 10-20  $\mu\text{m}$  increment step 2  $\mu\text{m}$  · 20-60  $\mu\text{m}$  increment step 5  $\mu\text{m}$  · Minimum setting value: 0.25  $\mu\text{m}$  Semi Automatic Rotary Microtome, Microtome Instrument Exporters, Microtome Instrument Manufacturers