

**Product Name :**  
Infant Electric Incubator**Product Code :**  
NMCBS15002**Description :**

Infant Electric Incubator We have established ourselves as a renowned manufacturer and supplier of Infant Electric Incubator. All our products are widely acclaimed by the customers for their corrosion & abrasion resistance, sturdy construction, low maintenance and easy handling. These products are designed and developed using premium raw material, which is sourced from the reliable vendors of the market. Other information: The controller is based on Micro computerized system Large bright displays are easy to monitor and for viewing from a distance Digital display system provided to measure the Infant temperature, Air temperature and to set the required temperature Feather touch keys provided for easy operation Bar graph LED system provided to indicate the heater output Baby probe / Air probe is made of highly sensitive sensor and interchangeable which does not need field calibration, the probe is made up of silicon sheathed cable, which allows easy cleaning. IV stand provided Storage drawers provided Large Baby tray provided Trendelenburg / fowler system provided Unit is powder coated Hygrometer provided to indicate the relative humidity in the baby compartment. Width: 1040 mm Depth: 600 mm Height: 1450 mm Weight: 45 kgs Approx Baby Bed Size: 550 x 390 mm Tilt: + / - 8o C Recommended Environment for operation: 24o C to 30o C. Humidity: 0 to 95 % R / H Electrical: 230 V AC + / - 10 % 50 Hz. Heater Wattage: 250 Watts Max. Over all Wattage: 450 Watts when heater, Fan is ON Fault indication and alarms for: Baby probe failure Air probe failure High temperature Low temperature High temperature cut off Heater failure Fan failure Power failure Infant Electric Incubator, Infant Incubator Exporters, Infant Incubator Manufacturers

---

Website: [SITE](#)

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**