

**Product Name :**  
Digital Ultrasonic Cleaner**Product Code :**  
NMCBS17003**Description :**

Digital Ultrasonic Cleaner Innovation is our forte and we strive to design and develop a precision engineered range of Digital Ultrasonic Cleaner. We design and develop all the medical lab products using high-grade raw material and latest technology to ensure that final output meets with international quality standards & norms. In addition to this, we provide customization as per the hospital and medical lab tenders specifications laid down by the customers.

**Technical Specification :**

Digital Ultrasonic Cleaner Specifications: Tank Dim.: 45cm x 15cm x 15cm Liter Cap.: 10 Ltrs U/S Power: 200 Watts U/S Freq.: 25 khz Input:180 Vac - 270 Vac Warranty: 1 Year S.M.P.S. ensures smooth working of cleaners within the voltage variation of 180-270 vac. Half the weight & size of the conventional machines Competitive rates efficient service backup of India Models can be made as per your specific requirements Heater ensures high efficiency in cleaning Cleaners are equipped with Imported Transducers for smooth and long life working. Industries Applicable:- Hospitals, Nursing Homes, Educational Institutes, Electronic Industry, Jewelry Industry, Hotel Industry, Pharmaceutical Industry & Automobile Industry Digital Ultrasonic Cleaner, Ultrasonic Cleaner Exporters, Ultrasonic Cleaner Manufacturers

