

**Product Name :**  
Refri - Freezer**Product Code :**  
NMCBS20001**Description :**

Refri Freezer Tender Suppliers Africa and Refri Freezer Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Refri Freezer Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Refri Freezer China, Refri Freezer Kenya.

Refri Freezer Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Refri Freezer Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point),

---

Suriname (Paramaribo)

**Technical Specification :**

Refri - Freezer Refri Freezers are most ideal for small blood banks/ laboratories where the upper refrigerator section can be used for storing blood bags vaccines, costly chemicals etc. as a Laboratory Refrigerator. · A TWO-IN-ONE/COMBO unit of a refrigerator and freezer. · Minimal foot- print required, a genuine space saving concept. · The high-efficiency compressor with internationally famous brand of Danfoss, accounts for faster and uniform cooling; · Unique refrigerating technique of refrigeration circulation is used for uniform cooling; · 100mm puff insulation effectively reduces loss of cold air; · German fan electromotor with internationally famous brand of Ebm is used. · Refri Freezer offers a wide variety of research and storage applications, such as low temperature scientific experiments, preservation of plasma, biomaterial, vaccine, biomedical product. · Most ideal for small blood banks/ laboratories where the upper refrigerator section can be used for storing blood bags vaccines, costly chemicals etc. as a Laboratory Refrigerator. · The lower LOW TEMPERATURE SECTION can be used to store plasma bags, testing of low temperature resistant properties of military products and has a wide applications for research institutes, electronic industries, chemical engineering industries, hospital, sanitation and antiepidemic stations, university laboratories, military industries. · Adjustable shelves in the refrigeration section are convenient to place articles; · Door with chromatic mirror inset, Alloy door knob with latent frame, Beautiful appearance, luxurious and elegant; · Foaming insulation heat plate inside and two-double sealed, prevent cold air leakage and promote excellent heat insulation. Refri - Freezer, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)