



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Platelet Agitator

Product Code :
NMCBS20008



Description :

Platelet Agitator Platelet agitator specially designed to maintain, platelets in continues motion at a fixed speed of 70 RPM · Maximum 36 standard Platelet bags can be accommodate at a time Platelet Agitator, Blood Bank Instrument Exporters, Blood Bank Instrument Manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,