

**Product Name :**  
CFL Phototherapy Unit**Product Code :**  
NMCBS26008**Description :**

CFL Phototherapy Unit CFL Provides the light intensity of 18 PW/cm<sup>2</sup>/nm (approx) at 45cm from Baby cot · Choice of 6 Blue or 4 Blue + 2 white CFL lamps for treatment and inspection · Height 120 cm to 160 cm from the base · Head Tilting Facility from -1800 + 1800 (Continuous) · 3 swivel casters fixed on sleek T-shaped base · Powder coated scratch and rust proof body · Special mirror Coated Reflector For better reflection CFL Phototherapy Unit, Infant Care Products Exporters, Infant Care Products Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)