



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Instrument Sterilizer

Product Code :
NMCBS26011



Description :

Instrument Sterilizer · A product of choice in the medical field · Seamless shell & lever operated Lid fitted with full proof mechanism control excessive steam escape and restricts condensate within the shell · Incredible synchronized maneuverability of lid, due to statistically perforated tray for flushing & entry of water. · Highly adoptable of your growing needs · Configuration design & material as per IS : 5022 · Our medical sterilizer, hospital sterilizer available in different sizes: 510 x 200 x 150 mm, 430 x 200 x 150 mm, 300 x 150 x 125 mm
Instrument Sterilizer, Infant Care Products Exporters, Infant Care Products Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,