

**Product Name :**  
Horizontal Sterilizer**Product Code :**  
NMALV10002**Description :**

Horizontal Sterilizer, Autoclave & Sterilizers Exporters, Autoclave & Sterilizers Manufacturers

**Technical Specification :**

The technique uses the method of gravity air removal to ensure uniform steam distribution and penetration

Steam trap equipped with air vent ensure maximum air removal. Useful for sterilization of glass containers i.e Ampoules, vials, glass bottles etc.

It is also useful for sterilization of glass wares, utensils, vessels, instruments, rubber articles linen etc

Hp the cycle 3 prester vacuum pulses-heat up-star hold- exhaust-post star-vacuum pulsing for drying.

This technique uses near 100% air removal from sterilization chamber v, which results in better temperature uniformity in the chamber because of the absence of air pockets i.e. Cold spots.

This technique increases the penetration of steam in porous products from which removal of air is difficult.

Leak test cycle for periodic testing of door integrity and pipeline connections & detecting leakage ampules

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**