

Product Name :
Stereo Acuity Test Kit

Product Code :
NMCBS36002



Description :

Stereo Acuity Test Kit
½ Designed to rapidly test for amblyopia and strabismus using gross thru fine Stereopsis (4800 to 20 seconds of arc)
½ Graded circle test from 400 seconds now down to 20 seconds
½ No monocular clues with new technology
½ New improved booklet 1005
½ Random Dot Stereopsis Test with LEA Symbols:
½ Designed to rapidly test for amblyopia and strabismus in early and non-readers and non-verbal children & adults
½ Expanded Random Dot LEA Symbols Test (500, 250, 125, 63 seconds of arc)
½ Graded circle test now down to 12.5 seconds with no monocular clues
½ New improved booklet 1010
½ Butterfly Stereopsis Test with LEA Symbols:
½ Designed to rapidly test for amblyopia and strabismus using gross through fine Stereopsis (2500 to 20 seconds of arc)
½ Graded circle test from 400 seconds now down to 20 seconds with no monocular clues
½ New improved booklet Stereo Acuity Test Kit, Medical Keratometer exporters, Medical Keratometer manufacturers