

**Product Name :**  
Multipara Patient Monitor**Product Code :**  
NMCBS38002**Description :**

Multipara Patient Monitor High bright LEDS display of NIBP, SPO2 waveform display Manual auto/continuous measurement of NIBP Adjustable audible and visual alarm Up to 600 groups NIBP data storage or SPO2 data up to 10 hours Suitable for adult, pediatric and neonatal patient SPO2: Measurement range: 0~100% Resolution: 1% Accuracy: 70~100% Accuracy: 70~100%+2%,0~69% unspecified Alarm range: 0~100% Range: 0~254bpm Resolution: 1bpm Dimension:200X230X75mm Weight is 3.5kg Power supply is AC 100~250V 50/60Hz Display: 7 segment LEDS for systolic, diastolic, mean pressures, SPO2 and pulse rate 2.4" color LCD waveform, bar-graphics, trend and system menu Multipara Patient Monitor, Multipara Monitor Systems exporters, Multipara Monitor Systems manufacturers