

**Product Name :**  
Digital Multipara Monitors**Product Code :**  
NMCBS38003**Description :**

Digital Multipara Monitors 8.4 = color TFT display Suitable for adult, pediatric and neonatal patient Basic parameters (ECG, SpO2, NIBP) in a durable case for bedside monitoring and transport monitoring Available and visual alarms with adjustable alarm ranges Net-workable data management and storage capacity ECG: Lead type: 5-lead Input: RA; LA; RL; LL; V Sweep Speed: 2.5mm/s, 25mm/s, 50mm/s Accuracy:  $\pm 1$ bpm or  $\pm 1\%$  which is greater Protection: withstand 4000ac/50hz voltage in isolation against electro surgical and defibrillation Have S-T detection and Arrhythmia analysis facility Have Alarm, audible and visual alarm, alarm events re-callable facility Standard Configuration: ECG NIBP SpO2 Optional: RESP2-TEMP 2-IBP ET(O2) Network Thermal Print Pacing Maker Digital Multipara Monitors, Multipara Monitor Systems exporters, Multipara Monitor Systems manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,