

Product Name :
M8B Multipara Monitor**Product Code :**
NMCBS38004**Description :**

M8B Multipara Monitor is a lightweight, compact, portable 10 inch wide color TFT LCD with Nellcor Oximax/EDAN SpO2, SpO2 pulse-tone modulation (Pitch Tone), multi-leads ECG waveforms (7 leads) display, arrhythmia analysis, multi-lead S-T segment, built-in rechargeable Lithium-ion battery, network capability, powerful data storage capacity (96 hours graphic and tabular trends of all parameters, 500 NIBP measurement, 60 alarm events). Optional accessories for adult, pediatrics and neonate.

Configuration: Model Standard Optional M8B 10 inch wide screen LCD, 5 ECG leads, RESP, NIBP, PR, EDAN SpO2, 2-TEMP, Li-ion Battery Thermal Recorder, 2-IBP, Nellcor OxiMax SpO2 M8B Multipara Monitor, Multipara Monitor Systems exporters, Multipara Monitor Systems manufacturers