

Product Name :
Multi Parameter Patient Monitor**Product Code :**
NMCBS38006**Description :**

Multi Parameter Patient Monitor
Multi Parameter Patient Monitor are High resolution 12.1" color TFT display
Multi Parameter Patient Monitor are Lightweight, compact and portable
Multi Parameter Patient Monitor are ECG, SPO2, NIBP, RESP, 2-TEMP PR
Multi Parameter Patient Monitor are Optional 2-IBP, ETCO2, and thermal printer
Multi Parameter Patient Monitor are Built-in rechargeable lithium battery
Multi Parameter Patient Monitor are ECG waveforms of 7-leads display on the same screen
Multi Parameter Patient Monitor are 72 hours graphic and tabular trends of all parameters
Multi Parameter Patient Monitor are 72 alarm events of all parameters recall
Multi Parameter Patient Monitor are 32 second full disclosure waveform review
Multi Parameter Patient Monitor are 500 NIBP measurement data can be storage and recall
Multi Parameter Patient Monitor are Data and waveforms color be adjustable
Multi Parameter Patient Monitor are Arrhythmia analysis and S-T segment analysis
Multi Parameter Patient Monitor are Suitable for adult, paediatric and neonatal patient
Multi Parameter Patient Monitor, Multipara Monitor Systems exporters, Multipara Monitor Systems manufacturers