

**Product Name :**  
Hand Held Single Channel ECG**Product Code :**  
NMCBS40001**Description :**

Hand Held Single Channel ECG Tender Suppliers Africa and Hand Held Single Channel ECG Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Hand Held Single Channel ECG Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Hand Held Single Channel ECG China, Hand Held Single Channel ECG Kenya.

Hand Held Single Channel ECG Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Hand Held Single Channel ECG Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward

---

Point), Suriname (Paramaribo)

**Technical Specification :**

1. Output system:High-resolution thermal-array(32 dots/mm),it needs not be adjusted.Recording frequency Response:150Hz(-3dB,+0.4dB).
2. The device can record exact single ECG waveform and remark. The mark includes:lead sign,ensitivity,paper speed,filter state.
3. Under automatic mode,just press the button once,it starts record procedure,which can enhance your work efficiency
4. The language of operation interface include (Chinese,English,French,Italian,Turkish,Spanish)
5. The keyboard's control is whole touch type,it's convenient to operate.
6. LCD can display the operation state, which is convenient and readable.
7. Classification:IEC class I,type CF
8. Power supply:Input AC: 100-240V, 50/60Hz Output DC: 9V DC:7.4V, Lithium rechargeable battery
9. This instrument can record 150 pieces of ECG waveform and 90 minutes continually under the best DC state
10. The figure of whole device is elegancy and gliding
11. According to defencence degree of deleterious fluid, this device is belong to common device.
12. According to the safe degree used under the condition with flammable anaesthesia gas mixed with air(or oxygen,nitrous oxide),this device belongs to the device which can't be used under the condition with flammable anaesthesia gas mixed with air(or oxygen,nitrous oxide).
13. Digital signal processor for effective inhibition of AC,EMG and DFT,receive the better ECG waveform.
14. According to the working mode class, this device belongs to continuous operation equipment.

ECG 80A Single channel Handheld ECG EKG machine with software electrocardiograph CE

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)