

**Product Name :**  
6 Channel ECG Machine**Product Code :**  
NMCBS40002**Description :**

6 Channel ECG Machine  
• Heart rate  
• Heart rhythm  
• For conduction abnormalities  
• Prior heart attack  
• Coronary artery disease  
• Abnormally thickening of heart muscle  
Features:  
• 6 Channel ECG Machine are 7" TFT, touch screen, led backlight  
• 6 Channel ECG Machine are Digital isolation technology and signal processing solution, digital filter  
• 6 Channel ECG Machine are Patient management, name and age and id.  
• 6 Channel ECG Machine are Detail analysis report.  
• 6 Channel ECG Machine are Multi-language options: Germany, France, Italian, Turkey, etc  
• 6 Channel ECG Machine are Built-in lithium rechargeable battery, wide thermal printing system  
• 6 Channel ECG Machine are USB and LAN socket.  
• 6 Channel ECG Machine are Memory: Built-in memory or mini SD Memory card. Store more than 1000 pieces of archive  
• 6 Channel ECG Machine are Auto-measurement, auto-analysis and auto-interpretation  
6 Channel ECG Machine, ECG Machines exporters, ECG Machines manufacturers