

**Product Name :**  
Physio Control Defibrillator**Product Code :**  
NMCBS41003**Description :**

Physio Control Defibrillator Owing to our experience, we have been successful in catering to the requirements of clients by offering quality Physio Control Defibrillator. All the products are widely acclaimed by the customers for their excellent performance, easy installation, low maintenance and longer service life. In addition to this, we provide customization as per the specifications laid down by the customers. Features: Physio Control Defibrillator are Simple, logical controls Physio Control Defibrillator are User-friendly displays Physio Control Defibrillator are Code summary Physio Control Defibrillator are Critical event record Physio Control Defibrillator are Space saving Physio Control Defibrillator are Vertical design Physio Control Defibrillator are Optional shock Physio Control Defibrillator are Advisory system Innovative features: · High and low energy · Select displays · Annotating recorder · Auxiliary paddles for “hands-free” DE fibrillation. · Designed for ease of maintenance · Diagnostics for simple test and calibration Physio Control Defibrillator, Electrical Defibrillator exporters, Electrical Defibrillator manufacturers