

Product Name :
Siemens 300 Ventilator**Product Code :**
NMCBS42002**Description :**

Siemens 300 Ventilator We are engaged in offering our clients with Siemens 300 Ventilator that is available in different models and sizes at market leading prices. These are designed in compliance with international quality norms & standards, using high-grade raw material and latest technology. In addition to this, we check all the products on certain parameters to ensure flawlessness. The Siemens 300A is comprised of 2 main components: Control Unit: Contains the electronic circuits necessary for control of ventilation (panel and knobs). Patient Unit: Continually measures flow and pressure with transducers. Features: Siemens 300 Ventilator are Auto mode Siemens 300 Ventilator are Neonate to adult patient range Siemens 300 Ventilator are Wide range of ventilation modes Siemens 300 Ventilator are Easy to use Siemens 300 Ventilator are O2 air gas module filters bacteria. Specifications: Ventilation modes: Controlled, Supported, Combined Patient range: Adult, Pediatric, Neonate Classification: IEC60601 - 1/EN60 601-1 Type B Dimensions: Patient unit: W 242 x D 370 x H 240 mm User interface: W 431 x D 150 x H 325mm Weight: Approx. 24kg Method of triggering: Flow and Pressure

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,