

**Product Name :**  
LP10 Ventilator**Product Code :**  
NMCBS42009**Description :**

LP10 Ventilator

- Pressure control knob for pressure-limited, time-cycled ventilation
- Decreased internal resistance reduces breathing effort in SIMV mode
- Continuous, low-maintenance operation with minimal supervision results in low cost of ownership
- Features and portability provide flexibility to meet a broad range of ventilation and patient requirements
- Internal and external battery operation provides convenience and security

Indications For Usage:

- A portable ventilator ideal for in hospital transport, sub acute facilities and long term, in-home care