

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Siemens 300A Ventilator

Product Code :
NMCBS42014



Description :

Siemens 300A Ventilator Auto mode · Neonate to adult patient range · Wide range of ventilation modes · O2 air gas module filters bacteria · Easy to use · Ventilation modes: Controlled, Supported, Combined · Patient range: Adult/Pediatric/Neonate · Classification: IEC60601-1/EN60 601-1 Type B · Dimensions: Patient unit: W 242 x D 370 x H 240mm · User interface: W 431 x D 150 x H 325mm · Method of triggering: Flow and Pressure Siemens 300A Ventilator, Medical Ventilator exporters, Medical Ventilator manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)