

Product Name :
Monochrome Display Finger Tip Oximeter**Product Code :**
NMCBS43002**Description :**

Monochrome Display Finger Tip Oximeter Monochrome Display Finger Tip Oximeter are Easy to use with one-key operation Monochrome Display Finger Tip Oximeter are High accuracy and repetition Monochrome Display Finger Tip Oximeter are Automatic power off without signal for 5 seconds Monochrome Display Finger Tip Oximeter are Display mode: LCD display Monochrome Display Finger Tip Oximeter are SPO2 display: 2-digit and 7-segment digital display Monochrome Display Finger Tip Oximeter are Pulse rate display: 3-digit and 7-segment digital display Monochrome Display Finger Tip Oximeter are Pulse wave intensity display: 6-segment bar-graph display Monochrome Display Finger Tip Oximeter are Low battery voltage indication: segment LCD signal

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,