

**Product Name :**  
Syringe Pump Dual Channel**Product Code :**  
NMCBS44001**Description :**

Syringe Pump Dual Channel Tender Suppliers Africa and Syringe Pump Dual Channel Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Syringe Pump Dual Channel Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Syringe Pump Dual Channel China, Syringe Pump Dual Channel Kenya.

Syringe Pump Dual Channel Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Syringe Pump Dual Channel Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines

---

(Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point),  
Suriname (Paramaribo)

**Technical Specification :**

Syringe Pump Dual Channel International standard syringe pump Compatible: Automatic syringe calibration, suitable for all brands Convenient: "Openable slider hook", convenient to clamp syringe plunger Powerful: Battery capable to support more than 6 hours for one channel continuous working Energy saving: Idle channel can be set to sleeping mode Function: Dual-channel with different flow rate Channel sleeping function 3-level occlusion alarm limits setting Limits setting Automatically reduce impact dose release after occlusion Delivered volume inquiry KVO flow rate Fast rate control "Purge" button activate safe protection Flow rate setting button lock Abundant optic and sound alarm AC and internal power indication Standard RS232 interface 500 history records Syringe Pump Dual Channel, Syringe Pump exporters, Syringe Pump manufacturers

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,