

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ceiling Mounted Operating Lights

Product Code :
NMCBS47012



Description :

Ceiling Mounted Operating Lights Ceiling Mounted Operating Lights are 19" dia metallic dome fitted with halogen bulb (1x 24V/150W) & heat absorbing colour correcting filters Ceiling Mounted Operating Lights are Provided with plastic diffusers for shadow reduction Ceiling Mounted Operating Lights are Spring balanced Ceiling Mounted Operating Lights are Provided with sterilizable handle and low voltage supply unit Ceiling Mounted Operating Lights are Light Intensity: 45,000 Lux + 10% Ceiling Mounted Operating Lights, X-Ray Machine & Instrument exporters, X-Ray Machine & Instrument manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,