

**Product Name :**  
Hot Plates

**Product Code :**  
NMALV10010



**Description :**

Hot Plates

**Technical Specification :**

Construction: Round Casted Plate.

Moc(material Of Construction): Body Made Of G.i. Duly Powder Coated.

Colour: Attractive Finished, Asthetic Look With Outer Body Duly

Coating: Powder Coated With White/ Siemens Grey And Panel Made Of Peacock Blue

Max Surface Temp.: 300 Deg C.

Heater: Embedded Nichrome Coil.

Test Method: Contain Very Wide Application From Drying, Heating To All Laboratories.

Type Of Hotplate: Electrically Heated Lab Type.