



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Pharmacy Lab Spectrophotometer

**Product Code :**  
NMCBS51005



**Description :**

Pharmacy Lab Spectrophotometer Wavelength Range 190-1100 nm Czerny Turner Spectral Bandwidth 1 nm Working Mode Basic/Quantative/Wavelength Scan/DNA Protein Test/Kinestics (T,A,C) Optical System Double beam (1200 lines/mm Grating) Scanning Speed Fast / Medium / Slow Photometric Range -0.3 - 3.5 A, 0 to 220 %T Wavelength Accuracy +0.3 nm Wavelength Reproducibility < 0.10 nm Photometric Accuracy + 0.3% T (0-100% T) +0.002 A (0-0.1 A) Photometric Reproducibility < 0.15 % T Stray Light < 0.1% T (Nal. 220nm, NaNO<sub>2</sub> 340 nm) Software CD 1 No. RS 232 Cable 1 No. Operation Manual 1 No. Detector Dual Silicon Photodiode Display 6 inches high light blue LCD Power AC:220V + 10%, 50 Hz. Dimensions 630 x 470 x 210 mm Weight 26 Kg Glass Cells 16 No. Quartz Cells 4 Nos. Fuse 3 Nos. Instrument Cover 1 No. Software CD 1 No. RS 232 Cable 1 No. Operation Manual 1 No. Pharmacy Lab Spectrophotometer, Double Beam UV-VIS Spectrophotometer exporters, Double Beam UV-VIS Spectrophotometer manufacturers



**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,