



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
UV Spectrophotometer

Product Code :
NMCBS52003



Description :

UV Spectrophotometer Single Beam Microprocessor UV/VIS Spectrophotometer Range 195 to 1000 nm with facility for automatic concentration, % Transmission, Absorbance & K Factor measurement. It has sample holder with 4 position selector control for 10 mm path length cuvettes. Facility has been provided for interface with any printer (RST 232C). Set of 2 quartz and 4 glass cuvettes is provided. SOFTWARE INCLUDED (Imported) UV Spectrophotometer, Single Beam UV-Vis Spectrophotometer exporters, Single Beam UV-Vis Spectrophotometer manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,