



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
UV Spectrophotometer

**Product Code :**  
NMCBS52003



**Description :**

UV Spectrophotometer Single Beam Microprocessor UV/VIS Spectrophotometer Range 195 to 1000 nm with facility for automatic concentration, % Transmission, Absorbance & K Factor measurement. It has sample holder with 4 position selector control for 10 mm path length cuvettes. Facility has been provided for interface with any printer (RST 232C). Set of 2 quartz and 4 glass cuvettes is provided. SOFTWARE INCLUDED (Imported) UV Spectrophotometer, Single Beam UV-Vis Spectrophotometer exporters, Single Beam UV-Vis Spectrophotometer manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)