

Product Name :
Digital Conductivity Meters**Product Code :**
NMCBS63001**Description :**

Digital Conductivity Meters Tender Suppliers Africa and Digital Conductivity Meters Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Digital Conductivity Meters Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Digital Conductivity Meters China, Digital Conductivity Meters Kenya.

Digital Conductivity Meters Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Digital Conductivity Meters Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines

(Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point),
Suriname (Paramaribo)

Technical Specification :

Digital Conductivity Meters Display 3 1/2 Digit LED 3 1/2 Digit LCD Range 5 Ranges 0 - 1000 Mohs/cm 2
Ranges 2 & 20 Mohs/cm Temperature Compensation Manual 0-50oC - Cell Constraint Adjustment Adjustable
Pre Adjustable Conductivity Cell Platinum DIP Platinum DIP Measuring Frequency 1000 Hz 1000 Hz Accuracy
+1% FSI +1% FSI Power 230v + 10% 1 x 9V Battery Dimensions 76 x 275 x 175 mm 85 x 28 x 165 mm Weight
2 kg. (Approx.) 0.8 kg. (Approx.) Accessories Conductivity Cell, Operation Manual, Dust Cover Digital
Conductivity Meters, Digital Conductivity TDS Meters exporters, Digital Conductivity TDS Meters manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,