

Product Name :
Microprocessor Based pH Meter**Product Code :**
NMCBS66001**Description :**

Microprocessor Based pH Meter Tender Suppliers Africa and Microprocessor Based pH Meter Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Microprocessor Based pH Meter Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Microprocessor Based pH Meter China, Microprocessor Based pH Meter Kenya.

Microprocessor Based pH Meter Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Microprocessor Based pH Meter Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the

Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point), Suriname (Paramaribo)

Technical Specification :

Microprocessor Based pH Meter Specifications Display 16 x 2 alpha numeric LCD display Range 0 to 14 pH 0 to + 1999.9 mV, Temp. 0 to 100oC Resolution .01 pH 0.1 mV Temp. 0.1 Deg. C Accuracy + 0.01 pH + 1 digit mV, Temp + 0.10 Calibration Auto / Manual Storage Upto 100 samples Power 230V + 50Hz Accessories pH Electrode, Temperature probe, Electrode Stand, Buffer solution, Dust cover, Manual

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [Medical Hospital Equipments Manufacturers,](#)