

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Plankton Net

Product Code :
NMCBS87010



Description :

Plankton Net Made of silk bolting net 875mm, net suspended by muslin reinforcement. There braided lead attached to swivel front. Rear of bag with 50mm opening for glass vial collection of plankton with 2 graduated collecting tubes 50ml. Capacity: BIO-1109.1 : Ring dia 250mm. BIO-1109.2 : Ring dia 300mm. BIO-1109.3 : Ring dia 375mm. Plankton Net, BIOLOGY exporters, BIOLOGY manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,