



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
AC Slip Ring Motor

Product Code :
NMCBS88007



Description :

AC Slip Ring Motor Our organization has gained recognition as the affluent manufacturer, trader, exporter and supplier of a high quality range of AC Slip Ring Motor. The offered motors are manufactured using high grade components and contemporary technology. We assure our clients that our range is quality tested on predefined parameters, in order to meet international standards. Apart from this, our prestigious clients can avail these motors as per their needs, at highly competitive price. Features: Reliable functioning Simple installation Sturdy construction AC Slip Ring Motor, Ac Motor exporters, Ac Motor manufacturers



Medical Hospital Equipments Manufacturers

Website: [SITE](#)