

**Product Name :**  
Baby Incubator**Product Code :**  
NMBCE10004**Description :**

Infant Incubator: Acrylic Single Wall Canopy with six port holes, Front Loading, Four Inlets for administration of O<sub>2</sub>, IV lines or Ventilator tubing. Sliding Baby Tray with External jerk free gear based Head Up / Down facility. X-Ray cassette sliding facility. Stainless Steel rust free Heater, blower sub assembly with changeable micro Air filters. Provision of an open shelf in the trolley. Humidity heater and Hygrometer. For Skin / Air / Manual Mode.

Control Modes: Skin / Air / Manual  
Range: Skin control: 34 to 38 C  
Air control: Room temp. to 39 C  
Resolution: 0.1 C  
Accuracy:  $\pm 0.2$  C  
Display: 1in bright skin & air temp.(7 segment LED)  
Status: On LCD  
Control Mode & Set temp  
Timer: Timer / APGAR Timer / Baby Attend  
Manual mode: Time Programmable  
Heater output settable Over temp. cut-off  
Alarm: Auto audio alarm mute, Mute time programmable  
Safety alarms: ON LCD  
The mentioned are default limit, can be set or alter as per requirement. Baby Incubator, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,