

Product Name :
Baby Incubator**Product Code :**
NMBCE10004**Description :**

Infant Incubator: Acrylic Single Wall Canopy with six port holes, Front Loading, Four Inlets for administration of O₂, IV lines or Ventilator tubing. Sliding Baby Tray with External jerk free gear based Head Up / Down facility. X-Ray cassette sliding facility. Stainless Steel rust free Heater, blower sub assembly with changeable micro Air filters. Provision of an open shelf in the trolley. Humidity heater and Hygrometer. For Skin / Air / Manual Mode.

Control Modes: Skin / Air / Manual
Range: Skin control: 34 to 38 Ce
Air control: Room temp. to 39 C
Resolution: 0.1 C
Accuracy: ± 0.2 C
Display: 1in bright skin & air temp.(7 segment LED)
Status: On LCD
Control Mode & Set temp
Timer: Timer / APGAR Timer / Baby Attend
Manual mode: Time Programmable
Heater output settable Over temp. cut-off
Alarm: Auto audio alarm mute, Mute time programmable
Safety alarms: ON LCD
The mentioned are default limit, can be set or alter as per requirement. Baby Incubator, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)