

Product Name :
Dental X-Ray Unit**Product Code :**
NMCBS101006**Description :**

Dental X-Ray Unit · The dental X-ray package is a true insurance of quality with all the benefits of the Advanced Technocracy Inc · An extremely high dielectric coefficient is obtained through the use of specially treated oil · The high voltage transformer comes with the best insulators · A special expansion chamber compensates for the dilation of the oil caused by heat from the tube · The protection against leakage radiation is stricter than the strictest international standards. The hooded tube with grid, compensates for voltage fluctuation · Input voltage: 230V AC single phase Timer: Electric with digital display of mAs Output: 8 mA - 70 KVP Stand: Floor model/Wall model-ray tube: 1.5mm focal Spot Cone: Parallel Dental X-Ray Unit, Dental Instrument exporters, Dental Instrument manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,