

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Curing Light

Product Code :
NMCBS101011



Description :

Curing Light · Power source: 230V~115V or 115V~(on request), 50/60Hz · Weight: 2.4 Kgs. · Dimension: W200 x D140 x H80(mm) Curing Light Function: · Complete curing restoration with a 75 watt, 12V, halogen lamp · Light guide with 360 degree rotation and can be autoclavable · 3 seconds interval to turn on or off the lamp controlled by computer to protect and lengthen use life of lamp · Individual curing time selection with beep tone signal · Automatical turning off fan after cool down period · Provides fuse and thermal cut-off switch to prevent overheating Curing Light, Dental Instrument exporters, Dental Instrument manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)