

Product Name :
Dental Double Walled Acrylizer**Product Code :**
NMCBS101014**Description :**

Dental Double Walled Acrylizer · Acrylizer Double walled, inside made of stainless steel and out side mild steel sheet painted with attractive stove enamel · Temperature is controlled by hydraulic thermostat from 5 c above ambient to 100 C with a sensitivity of + 5 C Acrylizer Consist Mechanical Timer can be set 1/2 Hour to 2 Hours, brass drain plug, control knob, On/Off switch · Acrylizer Inner Size : 300 X 250 X 350 mm Dental Double Walled Acrylizer, Dental Instrument exporters, Dental Instrument manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,