

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Automatic Dental Amalgamator

Product Code :
NMCBS101015



Description :

Automatic Dental Amalgamator · Automatic amalgamation with timer is a compact and functionally designed mixer which assures efficient and reliable trituration with a simple operation · Excellent durability and mixing power with 4 capsules · Dimensions: 13 cm. Width, 22cm. Length, 12 cm. Height, 3.5 kg Weight · Power 220/230 Volts AC Consumption: 45 W Automatic Dental Amalgamator, Dental Instrument exporters, Dental Instrument manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)