

Product Name :
Foetal Heart Monitor, Ultrasonic Doppler**Product Code :**
NMBCE10007**Description :**

Foetal Heart Monitor, Doppler: A high-performance Ultrasonic Doppler unit in shock-proof ABS moulded body for locating the position of Placenta or confirming the life and heartbeat of foetus.

Especially designed circuit design to perform well in a small size. › Enhanced performance - works even with obese patients or in case of hydramnios.

Low Battery & Mains indicators.

Small handheld size.

Light weight.

Rechargeable Batteries with built-in charger (Optional).

Operates on Mains also.

Earphone and recording sockets.

Technical Details:

Weight : 300 gm.

Height : 5.6in.

Width : 2.9in.

Depth : 1.4in.

Probe Transmitter Freq.: 2.25 MHz.

Battery Life : 300 X 1 minute.

Complies with : IEC 601-1.

Protection Class - I. Foetal Heart Monitor, Ultrasonic Doppler, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,