

Product Name :
Orthopedic Bed, Fowler**Product Code :**
NMCBS106004**Description :**

Orthopedic Bed, Fowler 1/2 Orthopaedic bed with mechanically operated back rest 0 to 70 1/2 & knee rest 0 to 35 1/2 all by Super Smooth crank mechanism. 1/2 Size: 2150 1/2 900 1/2 500 mm. 1/2 Epoxy Coated Steel frame work and 4 section perforated top. 1/2 Stainless Steel Orthopaedic attachment (without weight chain and weight hangers). 1/2 Head panel and foot panel equipped with safety lock and roller bumpers. 1/2 Mounted on protective stumps. 1/2 Provision for I.V. Rod. Freight saving knock-down construction Orthopedic Bed, Fowler, Hospital Bed exporters, Hospital Bed manufacturers